

VICAM – Vomitoxin Testing

User friendly lateral-flow
Don-V test delivers accurate
quantitative results as low as
0.20 ppm in a matter of minutes



Test the Vomitoxin in Your Fusarium Infected Grain

VICAM's Vertu® Lateral Flow Reader provides fast, easy and quantitative mycotoxin screening. Digital readings are clearly displayed on the screen of the Vertu Lateral Flow Reader, eliminating any guesswork about the results.

Vicam Don-V Test Strips with Vertu Reader

VICAM's Vertu® Lateral Flow Reader provides fast, easy and quantitative mycotoxin screening. Digital readings are clearly displayed on the screen of the Vertu Lateral Flow Reader, eliminating any guesswork about the results.

- Get results in less than 5 minutes (after filtration)
- No special training required
- Easily performed on-site or in the lab
- Accurate real-time data can be printed or downloaded to a computer for immediate use
- Wide Range – 0 to 10 ppm